HAMA

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

For further languages please visit: www.milenia-biotec.de or send an email to: info@milenia-biotec.de

Warnings and Precautions

- Read **these instructions for use** carefully before performing the test. Follow them step by step.
- Store reagents at 2 - 8° C.
- Do not use Test Unit if foil pouch is damaged.
- Protect Test Units from humidity.
- Do not interchange components of packages of different batch codes.
- All samples, used Test Units, pipettes, tips, and other contaminated material should be handled carefully and as potentially infectious.
- The disposal of waste materials must be carried out according to current local regulations.
- For *in vitro* diagnostic use only. These reagents are to be used only by certified medical laboratory personnel authorized by the laboratory.
- The reagents provided in these kits contain materials of human and/or animal origin. Whenever human blood product is required for the preparation of these reagents, approved methods are used to test the blood for the antibodies to HIV 1, HIV 2 and HCV, and for hepatitis B surface antigen, and results are found to be negative. However, no test method can offer complete assurance that infectious agents are absent. Therefore, users of reagents of these types must exercise extreme care in full compliance with safety precautions in the manipulation of these biological materials as if they were infectious.

Material Supplied, Storage, and Stability

REF: MQHM Z					
Component	Cat.-No.	Content	Preparation	Storage at	Shelf life
HAMA Test Unit; membrane coated with mouse IgG; mouse IgG conjugated to gold nanoparticles	MQSHM	5	ready- to -use	2 - 8 °C <i>Use immediately after opening the foil pouch!</i>	in the <u>closed</u> foil pouch until expiration date
Chase Buffer, Buffer, pH 7.2	MQBHM	1 vial à 0.6 mL (green)	ready- to -use	2 - 8 °C	until expiration date
MQHM Z: This bag contains 5 Evaluation Cards	MQAHMZ	5	ready- to -use		

Material Safety Data Sheets are available on request (look as well www.milenia-biotec.de).

Intended Use

For *in vitro* diagnostic use.

HAMA is designed for detection of human anti-mouse antibodies (IgG).

Materials Required But Not Supplied

Pipet and tips: 50 µL

Specimen Collection and Preparation

Serum.

Use tubes without any anticoagulant!

Avoid haemolysis; separate serum from blood cells within two hours. If there is no centrifuge available place serum sample tubes after blood collection (as described in sample tubes' IFU) vertically in a sample rack until the natural sedimentation has finished and a stable blood clot has been formed.

Volume Required

50 µL serum

Summary and Explanation

Monoclonal mouse-antibodies are used *in vivo* e.g. to prevent rejection of the graft or for immunotherapy in tumour patients, respectively. During the course of such therapy patient's can form antibodies against the mouse-antibodies (HAMAs). The formation of HAMAs can have an impact on the effectiveness of the therapy. In addition, HAMAs lead to troublesome cross-reactions in immunoassays, which can cause incorrect results.

Testprinciple

The HAMA is a lateral flow immunoassay designed for qualitative determination of human anti-mouse antibodies in human serum.

The sample will be pipetted in the sample application port of the test unit and the chase buffer is added immediately. The buffer forces the sample to migrate through the membrane of the test unit.

Human anti-mouse antibodies from patient's sample bind first to mouse antibodies which are conjugated to gold nanoparticles.

In a second step the human anti-mouse antibodies- from patient's sample- will be caught by mouse IgG antibodies, coated on the T-line of the nitrocellulose membrane.

In the presence of the HAMA-Conjugat-complexes a band becomes visible at the location of the test line (T). This test line- if it is positive- has to be the same or a darker shade of color than the given one on the provided evaluation card.

The surplus of gold particles continues to migrate through the membrane and is captured at the control line (C) by specific antibodies. A well visible line appears after incubation time.

Color's intensity of the test line is directly proportional to the HAMA concentration in the sample and increases during incubation time.

Procedure

Important notes:

- Before opening bring the test unit and sample buffer to room temperature (18 – 28 °C).
- Use the test unit immediately after opening. Do not store at room temperature!

	Steps	
1.	Sample 	Add 50 µL serum on the pad of test unit's sample application port.
2.	Buffer 	Add 1 drop chase buffer Note: Do not contaminate the vial with sample material! Avoid air bubbles!
3.	Incubation	10 minutes at room temperature (18 - 28°C).
4.	Evaluation: (visual) By the included Evaluation Card only. 	The result can only be interpreted visually by the provided Evaluation Card : The test is interpreted positive if the test line shows the same or a darker shade of color after incubation time than the given line in the evaluation card (s. results interpretation and functional control, p. 5) The result must be interpreted immediately after incubation time.

Results Interpretation

Put the colored line of the **Evaluation Card** (as shown below) **close beside** the test line of your test unit.

(See figures 1, 2)


Fig. 1


Fig. 2

The **result** will be evaluated **negative**, if the **Test Line-T** shows a **lighter-colored or weaker shade of color** than the given line on the Evaluation Card. If no **T-Line** is visible the result is also **negative**.

The **result** will be evaluated **positive**, if the **Test Line-T** shows the **same (Fig. 1) or a darker shade of color (Fig. 2)** than the given line on the Evaluation Card.

The **Control Line-C** must be visible. (See **Function Control** of the test unit.)


Note: Results should be verified concerning the entire clinical status of the patient. Additionally every decision for therapy should be taken individually.

Function Control


In any case, the Control line (C) has to appear!

It is only used as functional control and can not be used for the interpretation of the test line.
If the control line is not detectable after 10 minutes incubation time, the result is invalid and the analysis has to be repeated with a new test unit!

Valid results

Invalid results

Assay Characteristics

Sample Material:	50 μ L serum
Incubation:	10 minutes at room temperature
Sensitivity	97%
Specificity:	91%

Limitation / Interferences

Do not use plasma samples!

HAMA

Lateral Flow Immunoassay to detect human anti-mouse antibodies (HAMA)

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