Polycheck[®]

Scanner-based screening assay for the quantitative determination of allergen specific IgE

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

Materials Supplied, Storage and Stability

Components	Cat.-No.	Content	Preparation	Store at	Shelf Life
Polycheck® Pediatric Screen, Chip coated with 5 calibrators and 20 allergens: f13, f02, f01, f75, f35, f31, f03, f49, f14, f04, t03, g06, w06, d01, d02, e02, e01, e03, m03, m02	MCPPC 1(Z)	20 (10) Chips	ready to use	2 - 8 °C Protect from moisture! Store together with desiccant in a carefully sealed plastic bag	until the expiration date
Polycheck® Inhalation Screen, Chip coated with 5 calibrators and 20 allergens: t03, t02, t04, t07, g06, g12, w06, w09, d01, d02, e02, e01, e03, e06, e84, e82, m03, m02, m01, m06	MCIPC 1 (Z)	20 (10) Chips			
Polycheck® Food Screen, Chip coated with 5 calibrators and 20 allergens: f17, f13, f16, f20, f02, f01, f75, f78, f35, f85, f31, f25, f03, f24, f95, f49, f14, f04, f10, f05	MCFCPC 1 (Z)	20 (10) Chips			
Polycheck® Mediterranean Screen, Chip coated with 5 calibrators and 20 allergens: t03, t02, t04, t07, g06, g12, w01, w06, w09, w21, t09, d01, d02, e02, e01, e03, m03, m02, m01, m06	MCMPC 1 (Z)	20 (10) Chips			
Start Solution, buffered protein solution	MSPC	2 vials (1 vial) 5 mL, blue	ready to use	2 - 8 °C	30 days after opening, respectively until expiration date
anti-IgE Antibody, monoclonal (murine) antibody labeled with ligand	MAPC	2 vials (1 vial) 5 mL, green	ready to use	2 - 8 °C	30 days after opening, respectively until expiration date
Enzyme-Labeled anti-Ligand, Ligand conjugated to alkaline phosphatase	MEPC	2 vials (1 vial) 5 mL, white	ready to use	2 - 8 °C	30 days after opening, respectively until expiration date
Substrate Solution, 5'bromo-4'chloro-3' indolyl-phosphate-and 4' nitroblue tetrazolium, buffered	MBPS3	2 vials (1 vial) 5 mL, black	ready to use	2 - 8 °C protect from light!	30 days after opening, respectively until expiration date
Wash Buffer, phosphate buffer, pH 7.4	MWBPC	1 pouch	dissolve in 1 litre sterile demineralized water	2 - 8 °C	30 days after opening, respectively until expiration date

Material Safety Data Sheets are available on request (look as well www.milenia-biotec.de).

Materials Required

- Personal computer (Windows 98, Windows ME, Windows 2000, Windows XP), printer
- Flat- bed- scanner (requested types on demand)
- Biocheck Imaging Software (BIS) for patient-oriented analysis, calculation and report
- shaker (orbital shaker: approx. 120 rpm)
- variable pipette for 250 µl up to 1000 µL; Multipette (1 mL) to dispense the wash solution
- sterile demineralized water
- sterile-one-litre bottle (for diluted wash solution)

Specimen Collection and Preparation

Sera samples can be stored at 2-8 °C for up to 3 days. Storage at -20°C or lower is recommended for longer term storage. Avoid repeated freezing and thawing. Lipemic sera should be cleared by centrifugation.

Warnings and Precautions

All reagents of this kit are strictly intended for *in-vitro* diagnostic use only. Do not replace or mix up reagents from one and the same or several test packages. Reagents' bottles shelf-life (30 days) starts with opening the bottles. Opened bottles should be used beforehand opening a new bottle! Only the powdered wash buffer can be exchanged between several kit packages. The application should be used by staff particularly informed and educated in the practice of *in-vitro* diagnostics. Compliance with the prescribed protocol for test implementation is always recommended. Sample material of patients (e.g. serum or plasma) normally used in laboratory determinations are always considered potentially infectious. Samples of such high risk patient groups should always be marked clearly and if necessary be handled in safety work benches (e.g. lamina flow bench).

Method and Test Principle

The Polycheck[®] screening assay is an enzyme- immuno assay for the quantitative determination of allergen-specific IgE. There are 20 lines of relevant paediatric, food, or inhalant allergens coated together with 5 calibrators on a carrier membrane, which is located in the cavity of the Polycheck[®] chip. During incubation the patient's serum allergen-specific IgE binds to the corresponding allergens. Non-bound serum components are removed by washing. Ligand labeled anti-IgE antibodies bind to allergen-specific IgE. Unbound antibodies are removed by washing. Enzyme-labeled anti-ligand binds to the immune-complexes; surplus of enzyme conjugate is washed away. The substrate solution is added. Specifically bound enzymes convert the colourless substrate to a dark precipitate. The colour intensity of the (allergen-) lines is proportional to the respective allergen-specific IgE concentration in the patient's serum.

The Polycheck[®] chips can be evaluated using the Biocheck Imaging Software (BIS), a PC and a flat-bed scanner. The results are quantified in kU/L. This calculation is based on a number of calibrators, run on each chip.

Test Performance

Important notes:

- The required number of chips and test reagents should be at room temperature, before use.
 - Powdered wash buffer should be diluted with steril demineralized water at least 30 minutes prior to use. Use a steril bottle.
 - Do not allow the membrane of the test chips to dry during the assay is running.
 - All incubation steps are performed at room temperature (18-24°C) and with continual shaking.
1. Prepare a required number of Polycheck[®] allergen chips and mark them. Overlay allergen chips with **400 µL** Start Solution (blue) and incubate them for a few seconds. Tap the chips carefully on adsorbent paper.
 2. Add **400 µL** patient's serum and incubate for **45 minutes** on an orbital shaker.
 3. Decant sample, and wash **three times** with **1 ml** wash buffer.
For that place the chips onto the work bench and pipet carefully 1 ml wash buffer into the cavity of the chip. Decant washing solution from the wells into a sink. Afterwards, add two times more 1 ml wash buffer in the same way into the cavities. After this, tap the chips carefully on adsorbent paper again.
 4. Then, pipet **400 µL** wash buffer into every chip cavity and incubate for **5 minutes** on the shaker. Decant and wash again every chip **3x with 1 ml** wash buffer. Tap the chips carefully on adsorbent paper.
 5. Pipet **400 µL** anti-IgE Antibody (green) and incubate for **45 minutes** on the shaker.
 6. Decant and wash every chip **3x with 1 ml** wash buffer. Tap carefully chips on adsorbent paper.
 7. Add **400 µL** Enzyme-labeled Anti-Ligand (white) and incubate for **20 minutes** on the shaker.
 8. Decant and wash every chip **3x with 1 ml** wash buffer. Tap carefully on adsorbent paper.
 9. Pipet **400 µL** substrate solution (black) and incubate **20 minutes in the dark** on the shaker.
 10. Decant and wash every chip **3x with 1 ml** wash buffer. Tap carefully the chips on adsorbent paper.
 11. Let the Polycheck[®] Chips become **air-dry**. Before evaluation with a scanner and the BIS software the membranes **must be absolutely dry!**

The chips can be evaluated up to 3 month, if they will be stored at room temperature and protected from light and humidity.

Interpretation of Results

The evaluation of the Polycheck[®] allergen chips are operated with the assistance of a PC, the Biocheck Imaging Software (BIS), and a flat-bed scanner. The Software (BIS) analyzes, documents, and stores all patients' results in its database.

The Polycheck[®] Allergen chips are placed upside down onto the flat-bed-scanner's surface. Afterwards they are scanned and the data are processed by the software. - Please note orientation labelling on the surface.

The BIS-program verifies the technique of the chips, the validity of the calibrators, and ranges the intensity of greyscales and the position of the scanned band pattern to each allergen of the selected panel. BIS calculates allergen-specific IgE concentrations for each allergen on the membrane. A patient's report is generated for every sample. There are all allergens listed which have been tested in this panel with the present concentration of allergen specific IgE. The results were quantitatively reported in International Units [kU IgE/L]. The RAST-class score to the respective IgE concentrations, and its clinical interpretation is shown in the following table:

IgE [kU/L]	Class	Interpretation
< 0.35	0	No allergen specific antibody detectable
0.35 – 0.7	1	Very low antibody titer; often no clinical symptom- despite sensitization
0.7 – 3.5	2	Low antibody titer; sensitization, often clinical evidence in upper class ranges
3.5 – 17.5	3	Clear antibody titer; clinical evidence is mostly present
17.5 – 50	4	High antibody titer; always nearly with existing actuality clinical signs
50 – 100	5	Very high antibody titer
> 100	6	Extreme high antibody titer

Limitation of the Method

1. A definitive clinical diagnosis should not be based on the results of a single diagnostic method. It should only be made by a physician after all clinical and laboratory findings have been evaluated.
2. Allergenicity of food allergens may be modified by heating (cooking) and/or digestion. In this case the Polycheck[®] results may be inconsistent with the clinical findings.
3. Due to occurrence of common antigenic determinants in allergens of one natural substance family, the possibility of crossreactions cannot completely be excluded. With the Polycheck[®] system only the declared allergen panels can be covered, but not other (crossreacting) allergens!

Assay Characteristics

Sample material:	serum
Incubation time:	2 hours and 10 minutes at room temperature (18-24°C) on an orbital shaker
Lower detection limit:	0.15 kU/ L
Specificity	human IgE; no crossreaction with IgG, IgA, IgD, IgM
Calibration:	using single donor's and pool sera with defined allergen- specific IgE concentration
Method Comparison	skin tests and alternative solid-phase immunoassays

Short Instruction: Polycheck®

Steps	Solution	per Chip
Pipet	Start Solution (blue)	400 µl
Incubate a few seconds, decant and tap the chips on absorbant paper.		
Pipet	Sample	400 µl
Incubate 45 min at room temperature (RT) on an orbital shaker.		
Decant, wash 3x with 1mL wash buffer (per chip). Add 400 µL wash buffer and incubate 5 min at RT on the shaker. Decant, wash 3x with 1mL wash buffer, tap the chips on absorbant paper.		
Pipet	Wash Buffer	400 µl
Pipet	Ligand-labeled anti-IgE (green)	400 µl
Incubate 45 min at RT on an orbital shaker.		
Decant, wash 3x with 1mL wash buffer (per chip), tap the chips on absorbant paper.		
Pipet	Enzyme-labeled anti-Ligand (white)	400 µl
Incubate 20 min at RT on a shaker.		
Decant, wash 3x with 1mL wash buffer (per chip), tap the chips on absorbant paper.		
Pipet	Substrate Solution (black)	400 µl
Incubate 20 min at RT on a shaker in the dark		
Decant, wash 3x with 1mL wash buffer (per chip), tap the chips on absorbant paper.		
Air-dry the chips.		
Evaluate the dry chips by using Biocheck Imaging Software and a scanner.		

For detailed description of the procedure see page 4, also.

Available Panels:

Pediatric (MKPPC)		Inhalation (MKIPC)		Food (MKFPC)		Mediterranean (MKMPC)	
Peanut	f13	Birch Pollen	t03	Hazelnut	f17	Birch Pollen / <i>Betula alba</i>	t03
Milk	f02	Alder Pollen	t02	Peanut	f13	Alder Pollen / <i>Alnus spp.</i>	t02
Egg White	f01	Nutbush Pollen	t04	Walnut	f16	Nutbush Pollen / <i>Corylus avellana</i>	t04
Egg Yolk	f75	Oak Pollen	t07	Almond	f20	Oak Pollen / <i>Quercus robur</i>	t07
Potato	f35	Timothy Grass	g06	Milk	f02	Timothy Grass / <i>Phleum pratense</i>	g06
Carrot	f31	Cultivated Rye	g12	Egg White	f01	Cultivated Rye / <i>Secale cereale</i>	g12
Codfish	f03	Mugwort	w06	Egg Yolk	f75	Common Ragweed / <i>Ambrosia elatior</i>	w01
Apple	f49	English Plantain	w09	Casein	f78	Mugwort / <i>Artemisia vulgaris</i>	w06
Soybean	f14	<i>D. pteronyssinus</i>	d01	Potato	f35	English Plantain / <i>Plantago lanceolata</i>	w09
Wheat Flour	f04	<i>D. farinae</i>	d02	Celery	f85	Pellitory / <i>Parietaria judaica</i>	w21
Birch Pollen	t03	Dog Epithelia	e02	Carrot	f31	Olive Pollen / <i>Olea europea</i>	t09
Timothy Grass	g06	Cat Epithelia	e01	Tomato	f25	<i>Dermatophagoides pteronyssinus</i>	d01
Mugwort	w06	Horse Epithelia	e03	Codfish	f03	<i>Dermatophagoides farinae</i>	d02
<i>D. pteronyssinus</i>	d01	Guinea pig Epithelia	e06	Shrimp	f24	Dog Epithelia	e02
<i>D. farinae</i>	d02	Hamster Epithelia	e84	Peach	f95	Cat Epithelia	e01
Dog Epithelia	e02	Rabbit Epithelia	e82	Apple	f49	Horse Epithelia	e03
Cat Epithelia	e01	<i>A. fumigatus</i>	m03	Soybean	f14	<i>Aspergillus fumigatus</i>	m03
Horse Epithelia	e03	<i>C. herbarum</i>	m02	Wheat Flour	f04	<i>Cladosporium herbarum</i>	m02
<i>A. fumigatus</i>	m03	<i>P. notatum</i>	m01	Sesame	f10	<i>Penicillium notatum</i>	m01
<i>C. herbarum</i>	m02	<i>A. tenius</i>	m06	Cultivated Rye Flour	f05	<i>Alternaria tenuis</i>	m06

Polycheck®

Scanner-based screening assay for the quantitative determination of allergen specific IgE
English: Page 1-6

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