HIT

**Lateral Flow Immunoassay to detect
IgG antibodies against PF4/polyanion-complexes**
English: Page 1-6

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


Note: Significant changes are indicated by dotted lines in the margin.
In the end of the IFU you will find a table with causes of changes.

Warnings and Precautions

- Read **these instructions for use** carefully before performing the test. Follow them step by step.
- Store reagents at 2 - 8° C.
- Do not use Test Unit if foil pouch is damaged.
- Protect Test Units from humidity.
- Do not interchange components of packages of different batch codes.
- All samples, used Test Units, pipettes, tips, and other contaminated material should be handled carefully and as potentially infectious.
- The disposal of waste materials must be carried out according to current local regulations.
- For *in vitro* diagnostic use only. These reagents are to be used only by certified medical laboratory personnel authorized by the laboratory.
- The reagents provided in these kits contain materials of human and/or animal origin. Whenever human blood product is required for the preparation of these reagents, approved methods are used to test the blood for the antibodies to HIV 1, HIV 2 and HCV, and for hepatitis B surface antigen, and results are found to be negative. However, no test method can offer complete assurance that infectious agents are absent. Therefore, users of reagents of these types must exercise extreme care in full compliance with safety precautions in the manipulation of these biological materials as if they were infectious.

Materials Supplied, Storage and Stability

| REF: MQHIT 1 | | | | | |
|--|-----------------|-------------------------------|--------------------|--|---|
| Component | Cat.-No. | Content | Preparation | Storage at | Shelf life |
| HIT Test Unit; membrane coated with polyclonal (goat) anti-human IgG antibodies, anti-ligand antibodies conjugated to gold particles | MQSHIT | 20 | ready to use | 2 - 8 °C Use immediately after opening the foil pouch! | in the <u>closed</u> foil pouch until expiration date |
| Sample Buffer, Buffer, pH 7.4 | MQBHIT | 1 vial à 2 mL (white) | ready to use | 2 - 8 °C | until expiration date |
| MQHIT 1: This bag contains 20 Evaluation Cards | MQAHIT | 20 | | | |
| REF: MQHIT Z | | | | | |
| Component | Cat.-No. | Content | Preparation | Storage at | Shelf life |
| HIT Test Unit; membrane coated with polyclonal (goat) anti-human IgG antibodies, anti-ligand antibodies conjugated to gold particles | MQSHIT | 5 | ready to use | 2 - 8 °C Use immediately after opening the foil pouch! | in the <u>closed</u> foil pouch until expiration date |
| Sample Buffer, Buffer, pH 7.4 | MQBHITZ | 1 vial à 0.5 mL (white) | ready to use | 2 - 8 °C | until expiration date |
| MQHIT Z: This bag contains 5 Evaluation Cards | MQAHITZ | 5 | | | |

Material Safety Data Sheets are available on request (look as well www.milenia-biotec.de).

Intended Use

- For *in vitro* diagnostic use.
- Milenia® QuickLine HIT is designed for qualitative detection of IgG antibodies against PF4/polyanion-complexes in human citrated plasma **or** serum.

Materials Required But Not Supplied

Pipet and tips: 5 µL

Specimen Collection and Preparation

Serum or citrated plasma

Sample collection must be in conformity with the recommendations for haemostasis tests.

Blood (9 vol.) is collected in 0.109 M (i.e., 3.2 %) trisodium citrate anticoagulant (1 vol.) or in a tube without anticoagulant.

Centrifugation: Within two hours after the sample collection.

Do not use samples with visible fibrin strands. Such samples need to be centrifuged a second time until no fibrin strands are visible.

Volume Required

Citrated plasma or serum: 5 µL

Summary and Explanation

Heparin induced thrombocytopenia (HIT) is a life-threatening disease associated with exposure to unfractionated or (less commonly) low-molecular-weight heparin, drugs which are used very frequently for various indications. HIT II occurs in up to 5 % of patients on heparin. This prevalence varies with patient population and the nature of the disease and the type of heparin, thus the anticoagulant. Cardiac and major orthopedic surgeries expose a higher risk than other clinical situations when heparin is used.

HIT is caused by antibodies directed against a complex of platelet factor 4 (PF4) and polyanions such as heparin. The pathogenic and thrombogenic species is immunoglobulin IgG. There is evidence for a correlation between antibody concentration and clinical severity. Platelet activation and massively enhanced thrombin generation promote venous and/ or arterial thromboembolism. The rate of severe patient injury (e.g. amputations) and mortality is very high.

Patients with HIT II require alternative anticoagulation, but the drugs used in this setting have a high rate of hemorrhagic complications. There is no antidote available. Monitoring with standard methods is either not reliable or impossible. On top of this the drug costs are significantly higher than for heparins. Alternative anticoagulants should be reserved for patients with high clinical probability for HIT and with confirmed antibodies. Excluding HIT-antibodies is an important step towards selecting the most appropriate anticoagulation strategy.

Method and Test Principle

The Milenia® QuickLine HIT is a lateral flow immunoassay designed for qualitative determination of antibodies against PF4/polyanion-complexes in human citrated plasma or serum.

The sample is pipetted in the sample application port of the test unit and the sample buffer is added immediately. The buffer forces the sample to migrate through the membrane of the test unit.

Immobilized anti-human IgG on the nitrocellulose membrane binds patient's IgG-antibodies which are previously captured by the PF4/polyanion-complex which is detected by intensely colored gold nanoparticles.

The presence of PF4/polyanion-complexes becomes visible at a colored test line (T). This test line - if it is positive - has to be the same or a darker shade of color than the given one on the provided evaluation card.

The surplus of gold particles continues to migrate through the membrane and is captured at the control line (C) by specific antibodies. A well visible line appears after incubation time.

References

1. Greinacher A. Heparin-induced thrombocytopenia. J Thromb Haemost. 2009; 7 Suppl 1:9-12. Review.
2. Wilke T, Tesch S, Scholz A, Kohlmann T, Greinacher A. The costs of heparin-induced thrombocytopenia: a patient-based cost of illness analysis. J Thromb Haemost. 2009;7(5):766-73.
3. Bakchoul T, Giptner A, Najaoui A, Bein G, Santoso S, Sachs UJ. Prospective evaluation of PF4/heparin immunoassays for the diagnosis of heparin-induced thrombocytopenia. J Thromb Haemost. 2009;7(8):1260-5.
4. Sachs UJ, von Hesberg J, Santoso S, Bein G, Bakchoul T. Evaluation of a new nanoparticle-based lateral-flow immunoassay for the exclusion of heparin-induced thrombocytopenia (HIT). J Thromb Haemost. 2011; 106(6).
5. Warkentin TE, Greinacher A, Koster A, Lincoff AM. American College of Chest Physicians. Treatment and prevention of heparin-induced thrombocytopenia: American College of Chest Physicians Evidence-Based Clinical Practice Guidelines (8th Edition). Chest. 2008 Jun;133(6 Suppl):340S-380S.

Test Performance



Important notes:

- Do not interchange components of different lots!
- The test unit is sensitive to humidity. In case of damaged package do not use test unit!
- Before opening bring the test unit and sample buffer to room temperature (18 – 28 °C).
- Use the test unit immediately after opening. Do not store at room temperature!

| | Step | |
|----|--|--|
| 1. | Sample:  | 5 µL serum or citrated plasma on the pad of test unit's sample application port. |
| 2. | Buffer | Add 2 drops sample buffer Note: Do not contaminate the vial with sample material! Avoid air bubbles! |
| 3. | Incubation | 10 minutes at room temperature (18 - 28°C). |
| 4. | Evaluation: (visual) By the included Evaluation Card only.  | The result can only be interpreted visually by the provided Evaluation Card : The test is interpreted positive if the test line shows the same or a darker shade of color after incubation time than the given line in the evaluation card (s. results interpretation and functional control) The result must be interpreted immediately after incubation time. |

Results Interpretation

Put the colored line of the **Evaluation Card** (as shown below) **close beside** the test line of your test unit. (See figures 1, 2)

Figures 1

Figures 2

The **result** will be evaluated **negative**, if the **Test Line-T** shows a **lighter-colored or weaker shade of color** than the given line on the Evaluation Card. If no **T-Line** is visible the result is also **negative**.

The **result** will be evaluated **positive**, if the **Test Line-T** shows the **same (Fig. 1) or a darker shade of color (Fig. 2)** than the given line on the Evaluation Card.

The **Control Line-C** must be visible. (See **Function Control** of the test unit.)


Note: Results should be verified concerning the entire clinical status of the patient. Additionally every decision for therapy should be taken individually.

Function Control


In any case, the Control line (C) has to appear!

It is only used as functional control and can not be used for the interpretation of the test line.
If the control line is not detectable after 10 minutes incubation time, the result is invalid and the analysis has to be repeated with a new test unit!

Valid results

Invalid results

Assay Characteristics

| | |
|------------------|---|
| Sample material: | 5 µL serum or citrated plasma |
| Incubation time: | 10 minutes at room temperature (18 – 28 °C) |
| Sensitivity: | 100% |
| Specificity: | 93% |

Limitation / Interferences

Lipemia, bilirubin, and haemolysis do not show any interference influence on the results.

Lateral Flow Immunoassay to detect
IgG antibodies against PF4/polyanion-complexes
English: Page 1-6

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Referenzen

1. Greinacher A. Heparin-induced thrombocytopenia. J Thromb Haemost. 2009; 7 Suppl 1:9-12. Review.
2. Wilke T, Tesch S, Scholz A, Kohlmann T, Greinacher A. The costs of heparin-induced thrombocytopenia: a patient-based cost of illness analysis. J Thromb Haemost. 2009;7(5):766-73
3. Bakchoul T, Giptner A, Najaoui A, Bein G, Santoso S, Sachs UJ. Prospective evaluation of PF4/heparin immunoassays for the diagnosis of heparin-induced thrombocytopenia. J Thromb Haemost. 2009;7 (8):1260-5.
4. Sachs UJ, von Hesberg J, Santoso S, Bein G, Bakchoul T. Evaluation of a new nanoparticle-based lateral-flow immunoassay for the exclusion of heparin-induced thrombocytopenia (HIT). J Thromb Haemost. 2011; 106(6).
5. Warkentin TE, Greinacher A, Koster A, Lincoff AM; American College of Chest Physicians. Treatment and prevention of heparin-induced thrombocytopenia: American College of Chest Physicians Evidence-Based Clinical Practice Guidelines (8th Edition). Chest. 2008 Jun;133(6 Suppl):340S-380S.

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