## Total IgE

Lateral Flow Immunoassay for the quantitative evaluation of human Total IgE (TIgE); Interpretation with the Milenia®POCScan Reader  
English: Page 1-11

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**Note:**


Effective from lot 013 and higher:

**Transfer pipettes**  
**Lancets**  
**Cleansing swabs**  
**Adhesive tapes**

are no longer components of this kit!

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

**Note** : Significant changes are indicated by dotted lines in the margin.  
In the end of the IFU you will find a table with causes of changes.

## Warnings and Precautions

- Read **these instructions for use** carefully before performing the test. Follow them step by step.
- Store reagents at 2 - 8° C.
- Do not use Test Unit if foil pouch is damaged.
- Protect Test Units from humidity.
- Do not interchange components of packages of different batch codes.
- All samples, used Test Units, pipettes, tips, and other contaminated material should be handled carefully and as potentially infectious.
- The disposal of waste materials must be carried out according to current local regulations.
- For *in vitro* diagnostic use only. These reagents are to be used only by certified medical laboratory personnel authorized by the laboratory.
- The reagents provided in these kits contain materials of human and/or animal origin. Whenever human blood product is required for the preparation of these reagents, approved methods are used to test the blood for the antibodies to HIV 1, HIV 2 and HCV, and for hepatitis B surface antigen, and results are found to be negative. However, no test method can offer complete assurance that infectious agents are absent. Therefore, users of reagents of these types must exercise extreme care in full compliance with safety precautions in the manipulation of these biological materials as if they were infectious.

## Materials Supplied, Storage and Stability

Component	Cat.-No.	Content	Preparation	Storage at	Shelf life
<b>Total IgE Test Unit,</b> membrane coated with anti-ligand and monoclonal (mouse) anti-hIgE antibodies conjugated to gold particles	MQSTE	20	ready to use	2 - 8 °C Use <b>immediately</b> after opening the foil bag!	in the closed foil pouch until expiration date
<b>Assay Buffer,</b> buffer containing ligand labeled polyclonal (goat) anti-IgE antibodies	MQTEB	1 vial à 1.5 mL	ready to use	2 - 8 °C	until expiration date
<b>Certificate of Analysis</b>	CoA-MQTE	1			

Material Safety Data Sheets are available on request (look as well [www.milenia-biotec.de](http://www.milenia-biotec.de)).

## Intended Use

For *in vitro* diagnostic use.

Milenia® QuickLine Total IgE is designed for evaluation of human Total IgE in whole blood.

For quantitative evaluation use Milenia® POCScan Reader.

## Materials Requires But Not Supplied

**Note: Transfer pipettes, lancets, cleansing swabs and adhesive tapes are no longer components of this kit!**

Pipet and tips: 20µL

## Specimen Collection and Preparation

**Note: Transfer pipettes, lancets, cleansing swabs and adhesive tapes are no longer components of this kit!**

The Milenia® QuickLine Total IgE test runs with 20 µL capillary blood.

## Volume Required

20 µL capillary blood

## Method and Test Principle

Milenia<sup>®</sup> QuickLine Total IgE is a lateral flow immunoassay. It is designed to detect total immunoglobulin E (IgE) from whole blood. With Milenia<sup>®</sup> POCScan Reader

Immediately, after pipetting the sample (whole blood) into test unit's sample application port, add one drop of assay buffer.

IgE antibodies from patient's sample will be bound specifically to ligand-labeled polyclonal anti-IgE antibodies present in the assay buffer. These immune complexes migrate into the test strip and associate there with the gold conjugate, which are coated with monoclonal anti-IgE antibodies. The gold particles overflow the test line (T) where an anti-ligand is immobilized. Only the immune complexes consisting of ligand labeled anti-IgE antibodies, patient's IgE and the gold conjugate will bind to the anti-ligand. Over time a reddish test line is generated. Gold conjugate-specific antibodies, which are printed as a line on the membrane (control line, C), catch unbound gold conjugate; so a clearly visible line is generated during the test time.

Interpretation of results is only possible with the Milenia<sup>®</sup> POCScan Reader, which allows an objective interpretation and documentation of the result.

The total IgE concentration correlates directly with the color intensity of the test line (T). The Milenia<sup>®</sup> POCScan Reader measures the color reflection of the test line and calculates the result according a stored lot specific standard curve.

# Milenia® POCScan Reader Overview




**Battery position and replacement:**









- Please use Nickel-Metal-Hybride (NiMH) batteries only, which offer longer durability.
- **ATTENTION: Recharging of batteries is working only, if POCScan reader is switched ON!**

## Test Performance

### Important notes:


- The test is sensitive to humidity. **In the case of damaged package do not use test unit!**
- Before opening bring test unit and assay buffer to room temperature (18 – 28 °C).  
Use the test unit immediately after opening. Do not use kit reagents from different lots!
- **Please read the instruction carefully prior performing the test; follow the procedure step by step and notice the sequence!**
- **Start measurement immediately after adding the assay buffer to the patient sample in the sample port!** The incubation time of 7 minutes (420s) is controlled by the instrument.
- For identification the operating elements of Milenia® POCScan Reader see “Overview” on page 5.




A) PREPARATIONS - 1		
Step		* Remark: All figures (fig.) are example displays only!
1.		Press “ON/OFF” button  for at least two seconds To turn on the device. A self-test runs automatically. Menu „SCAN“ is automatically selected.
2.		Take the test unit out of the aluminium foil bag.  Open the drawer of Milenia® POCScan Reader and place the test unit into the drawer.
3.		Start “SCAN” menu. For that press “ENTER” button.
4.		Enter “USER ID” and “SAMPLE ID” with respective buttons: (UP/DOWN/BACK/FORWARD/ENTER).  Select “Next” and press  .

	Step	A) PREPARATIONS - 2
5.		<p>Compare lot number and product code displayed by the Milenia® POCScan Reader, with the lot number and the product code of the outer box label of the kit.</p> <p>If they are identical, select “Yes” and press </p> <p>The display automatically switches to <b>step 9</b>.</p>
6.		<p>If the displayed lot number is <u>not</u> identical with the lot number of the used kit, lot number and respective curve parameters need to be changed.</p> <p>In this case select “No” and press </p> <p>The Menu “Enter coefficients” (step 7a) opens.</p>
7a.		<p>Select Lot, Prm A, Prm F and Prm G with respective buttons  and change all respective data according the values listed in the “<b>Certificate of Analysis</b>” of the currently used Milenia® QuickLine Total IgE kit.</p> <p>Select “Next” and press </p>
7b.		<p>Select Prm H, Prm J, and Prm K with respective buttons </p> <p>and change all respective data according the values listed in the “<b>Certificate of Analysis</b>” of the currently used Milenia® QuickLine Total IgE kit.</p> <p>Select “Next” and press  The display switches back to the Scan Mode</p>
7c.		<p>Compare displayed Checksum with the Checksum listed in the “<b>Certificate of Analysis</b>” of the currently used Milenia® QuickLine Total IgE kit.</p> <p>If the checksums are identical, select “Yes” and press </p>
8.		<p>If Checksum is <u>not</u> identical, please select “No”, press  and accurately check in <b>step 7a-c</b> the different parameters of the “<b>Certificate of Analysis</b>”.</p>
9.		<p>The menu “START” appears.</p> <p><b>Do <u>not</u> activate this menu yet!</b></p> <p><b>Now, draw blood from patient’s fingertip. And continue with part “TEST (b)”</b></p>

Now, draw blood from patient's fingertip and continue with TEST (B) as written below:

B) TEST		
	Step	
1.		Pipet 20µL of patient's capillary blood into the application port of the test unit. (See fig.).
2.		Add <u>immediately</u> 1 drop assay buffer. (See fig.).  Note: Do not contaminate the vial with sample material!  Avoid air bubbles!
3.		Close the drawer of the Milenia® POCScan Reader!  START  measurement  Continue with part "C, MEASUREMENT".

<b>C) MEASUREMENT</b>		
After starting measurement.....		
Step		
1.		<p>Press  Enter button immediately to start Scan Mode.</p>
2.		<p>Incubation time of 420 seconds (7 min) is counting down (Test strip delay time).</p> <p>(Scan runs automatically).</p>
3.		<p><b>After 7 minutes the result of Milenia® QuickLine Total IgE will be displayed.</b>  <b>The result must not be interpreted later!</b></p> <p>“Next test” is selected by; pressing  to return to “Main Menu”.</p> <p>Now, the system is ready for a new measurement, <u>or</u> switch off the system by pressing the “ON/OFF” button  for two seconds.</p>

## Interpretation of Results

The Milenia® QuickLine Total IgE can be evaluated only with the Milenia® POCScan Reader. The calculated value is expressed in U/mL.

## Internal Quality Control

Only valid results are shown in the display of Milenia® POCScan Reader. If the result is shown up as invalid, the analysis has to be repeated with a new test unit! For further information see: “Milenia® POCScan Reader Manual.”

## Assay Characteristics

Sample material:	20 µL capillary blood
Incubation time:	7 minutes (420s) at room temperature (18 – 28 °C)
Sensitivity:	15 U/mL
Dynamic range:	30 U/mL to 2,000 U/mL
Specificity:	Highly specific reaction with human immunoglobulin E
Method Comparison:	High correlation with a world wide established laboratory method for determination of Total IgE from serum/plasma
Precision:	Interassay-CVs ≤ 15% over entire dynamic range
Interferences:	Lipemia, bilirubin and haemolysis have no interference influence on the results

**Total IgE**

**Lateral Flow Immunoassay for the quantitative evaluation of human Total IgE (TIgE); Interpretation with the Milenia®POCScan Reader**  
English: Page 1-11

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